

Shamir Attitude III® SV

Product Availability and Markings

Material	Index	Power Range	Nominal Base Curve*	Color	Invisible Marking
Hard Resin	1.498	-6.50 to +6.00	4, 5, 6, 7, 8	Clear	15 Ⓢ AT ₃ 250
Transitions®	1.5	-6.50 to +6.00	4, 5, 6, 7, 8	Gray & Brown	TR Ⓢ AT ₃ 250
Transitions® XTRActive™	1.5	-6.50 to +6.00	4, 5, 6, 7, 8	Gray	X5 Ⓢ AT ₃ 250
Polarized	1.498	-6.50 to +6.00	4, 5, 6, 7, 8	Gray & Brown	PZ Ⓢ AT ₃ 250
Drivewear®	1.5	-6.50 to +6.00	4, 5, 6, 7, 8	Drivewear®	DW Ⓢ AT ₃ 250
DLC™ (Trivex)	1.53	-6.50 to +6.00	4, 5, 6, 7, 8, 9, 10	Clear	TX Ⓢ AT ₃ 250
DLC™ (Trivex) Transitions®	1.53	-6.50 to +6.00	4, 5, 6, 7, 8, 9, 10	Gray & Brown	TTX Ⓢ AT ₃ 250
DLC™ (Trivex) Transitions® XTRActive™	1.53	-6.50 to +6.00	4, 5, 6, 7, 8, 9, 10	Gray	TXX Ⓢ AT ₃ 250
DLC™ (Trivex) Polarized	1.53	-6.50 to +6.00	4, 5, 6, 7, 8, 9, 10	Gray & Brown	PTX Ⓢ AT ₃ 250
1.56	1.56	-6.50 to +6.00	4, 5, 6, 7, 8, 9, 10	Clear	56 Ⓢ AT ₃ 250
Polycarbonate	1.586	-6.50 to +6.00	4, 5, 6, 7, 8, 9, 10	Clear	PC Ⓢ AT ₃ 250
Polycarbonate Transitions®	1.586	-6.50 to +6.00	4, 5, 6, 7, 8, 9, 10	Gray & Brown	PT Ⓢ AT ₃ 250

* The true base curve can be selected from a range of curves that has been defined in the Freeform Software.



ReCreating Perfect Vision®

Shamir Attitude III® SV

Product Availability and Markings

Material	Index	Power Range	Nominal Base Curve*	Color	Invisible Marking
Polycarbonate Transitions® XTRActive™	1.586	-6.50 to +6.00	4, 5, 6, 7, 8, 9, 10	Gray	PX Ⓢ AT ₃ 250
Polycarbonate Polarized	1.586	-6.50 to +6.00	4, 5, 6, 7, 8, 9, 10	Gray & Brown	PPZ Ⓢ AT ₃ 250
Polycarbonate Drivewear®	1.586	-6.50 to +6.00	4, 5, 6, 7, 8, 9, 10	Drivewear®	PW Ⓢ AT ₃ 250
Superlite™ 1.6	1.6	-6.50 to +6.00	4, 5, 6, 7, 8, 9, 10	Clear	16 Ⓢ AT ₃ 250
Superlite 1.6 Tribrid™	1.6	-6.50 to +6.00	4, 5, 6, 7, 8, 9, 10	Clear	TB Ⓢ AT ₃ 250
Superlite 1.6 Transitions®	1.6	-6.50 to +6.00	4, 5, 6, 7, 8, 9, 10	Gray & Brown	T6 Ⓢ AT ₃ 250
Superlite 1.6 Transitions® XTRActive™	1.6	-6.50 to +6.00	4, 5, 6, 7, 8, 9, 10	Gray	X6 Ⓢ AT ₃ 250
Superlite 1.6 Polarized	1.6	-6.50 to +6.00	4, 5, 6, 7, 8, 9, 10	Gray & Brown	PZ6 Ⓢ AT ₃ 250
Superlite 1.67	1.67	-6.50 to +6.00	4, 5, 6, 7, 8, 9, 10	Clear	67 Ⓢ AT ₃ 250
Superlite 1.67 Transitions®	1.67	-6.50 to +6.00	4, 5, 6, 7, 8, 9, 10	Gray & Brown	T67 Ⓢ AT ₃ 250
Superlite 1.67 Transitions® XTRActive™	1.67	-6.50 to +6.00	4, 5, 6, 7, 8, 9, 10	Gray	X67 Ⓢ AT ₃ 250
Superlite 1.67 Polarized	1.67	-6.50 to +6.00	4, 5, 6, 7, 8, 9, 10	Gray & Brown	PZ67 Ⓢ AT ₃ 250

* The true base curve can be selected from a range of curves that has been defined in the Freeform Software.



SHAMIR

ReCreating Perfect Vision®