

Shamir Attitude III® Sport

Product Availability and Markings

| Material | Index | Power Range | Nominal Base Curve* | Add | Color | Fitting Height | Invisible Marking |
|--|-------|----------------|----------------------|--------------|--------------|----------------|----------------------|
| Hard Resin | 1.498 | -6.50 to +5.25 | 4, 5, 6, 7, 8 | 0.75 to 3.50 | Clear | 18mm | 15 ☉ ☽ 250 |
| Transitions® | 1.5 | -6.50 to +5.25 | 4, 5, 6, 7, 8 | 0.75 to 3.50 | Gray & Brown | 18mm | TR ☉ ☽ 250 |
| Transitions® XTRActive™ | 1.5 | -6.50 to +5.25 | 4, 5, 6, 7, 8 | 0.75 to 3.50 | Gray | 18mm | X5 ☉ ☽ 250 |
| Polarized | 1.498 | -6.50 to +5.25 | 4, 5, 6, 7, 8 | 0.75 to 3.50 | Gray & Brown | 18mm | PZ ☉ ☽ 250 |
| Drivewear® | 1.5 | -6.50 to +5.25 | 4, 5, 6, 7, 8 | 0.75 to 3.50 | Drivewear® | 18mm | DW ☉ ☽ 250 |
| DLC™ (Trivex) | 1.53 | -6.50 to +6.00 | 4, 5, 6, 7, 8, 9, 10 | 0.75 to 3.50 | Clear | 18mm | TX ☉ ☽ 250 |
| DLC™ (Trivex) Transitions® | 1.53 | -6.50 to +6.00 | 4, 5, 6, 7, 8, 9, 10 | 0.75 to 3.50 | Gray & Brown | 18mm | TTX ☉ ☽ 250 |
| DLC™ (Trivex) Transitions® XTRActive™ | 1.53 | -6.50 to +6.00 | 4, 5, 6, 7, 8, 9, 10 | 0.75 to 3.50 | Gray | 18mm | TXX ☉ ☽ 250 |
| DLC™ (Trivex) Polarized | 1.53 | -6.50 to +6.00 | 4, 5, 6, 7, 8, 9, 10 | 0.75 to 3.50 | Gray & Brown | 18mm | PTX ☉ ☽ 250 |
| 1.56 | 1.56 | -6.50 to +6.00 | 4, 5, 6, 7, 8, 9, 10 | 0.75 to 3.50 | Clear | 18mm | 56 ☉ ☽ 250 |
| Polycarbonate | 1.586 | -6.50 to +6.00 | 4, 5, 6, 7, 8, 9, 10 | 0.75 to 3.50 | Clear | 18mm | PC ☉ ☽ 250 |
| Polycarbonate Transitions® | 1.586 | -6.50 to +6.00 | 4, 5, 6, 7, 8, 9, 10 | 0.75 to 3.50 | Gray & Brown | 18mm | PT ☉ ☽ 250 |

* The true base curve can be selected from a range of curves that has been defined in the Freeform Software.



ReCreating Perfect Vision®

Shamir Attitude III® Sport

Product Availability and Markings

| Material | Index | Power Range | Nominal Base Curve* | Add | Color | Fitting Height | Invisible Marking |
|--|-------|----------------|----------------------|--------------|--------------|----------------|--|
| Polycarbonate Transitions® XTRActive™ | 1.586 | -6.50 to +6.00 | 4, 5, 6, 7, 8, 9, 10 | 0.75 to 3.50 | Gray | 18mm | PX   250 |
| Polycarbonate Polarized | 1.586 | -6.50 to +6.00 | 4, 5, 6, 7, 8, 9, 10 | 0.75 to 3.50 | Gray & Brown | 18mm | PPZ   250 |
| Polycarbonate Drivewear® | 1.586 | -6.50 to +6.00 | 4, 5, 6, 7, 8, 9, 10 | 0.75 to 3.50 | Drivewear® | 18mm | PW   250 |
| Superlite™ 1.6 | 1.6 | -6.50 to +6.00 | 4, 5, 6, 7, 8, 9, 10 | 0.75 to 3.50 | Clear | 18mm | 16   250 |
| Superlite 1.6 Tribrid™ | 1.6 | -6.50 to +6.00 | 4, 5, 6, 7, 8, 9, 10 | 0.75 to 3.50 | Clear | 18mm | TB   250 |
| Superlite 1.6 Transitions® | 1.6 | -6.50 to +6.00 | 4, 5, 6, 7, 8, 9, 10 | 0.75 to 3.50 | Gray & Brown | 18mm | T6   250 |
| Superlite 1.6 Transitions® XTRActive™ | 1.6 | -6.50 to +6.00 | 4, 5, 6, 7, 8, 9, 10 | 0.75 to 3.50 | Gray | 18mm | X6   250 |
| Superlite 1.6 Polarized | 1.6 | -6.50 to +6.00 | 4, 5, 6, 7, 8, 9, 10 | 0.75 to 3.50 | Gray & Brown | 18mm | PZ6   250 |
| Superlite 1.67 | 1.67 | -6.50 to +6.00 | 4, 5, 6, 7, 8, 9, 10 | 0.75 to 3.50 | Clear | 18mm | 67   250 |
| Superlite 1.67 Transitions® | 1.67 | -6.50 to +6.00 | 4, 5, 6, 7, 8, 9, 10 | 0.75 to 3.50 | Gray & Brown | 18mm | T67   250 |
| Superlite 1.67 Transitions® XTRActive™ | 1.67 | -6.50 to +6.00 | 4, 5, 6, 7, 8, 9, 10 | 0.75 to 3.50 | Gray | 18mm | X67   250 |
| Superlite 1.67 Polarized | 1.67 | -6.50 to +6.00 | 4, 5, 6, 7, 8, 9, 10 | 0.75 to 3.50 | Gray & Brown | 18mm | PZ67   250 |

* The true base curve can be selected from a range of curves that has been defined in the Freeform Software.



ReCreating Perfect Vision®

Shamir Attitude III® Fashion

Product Availability and Markings

| Material | Index | Power Range | Nominal Base Curve* | Add | Color | Fitting Height | Invisible Marking |
|---------------------------------------|-------|----------------|----------------------|--------------|--------------|----------------|----------------------|
| Hard Resin | 1.498 | -6.50 to +5.25 | 4, 5, 6, 7, 8 | 0.75 to 3.50 | Clear | 15mm and 18mm | 15 ☉ ☹ 250 |
| Transitions® | 1.5 | -6.50 to +5.25 | 4, 5, 6, 7, 8 | 0.75 to 3.50 | Gray & Brown | 15mm and 18mm | TR ☉ ☹ 250 |
| Transitions® XTRActive™ | 1.5 | -6.50 to +5.25 | 4, 5, 6, 7, 8 | 0.75 to 3.50 | Gray | 15mm and 18mm | X5 ☉ ☹ 250 |
| Polarized | 1.498 | -6.50 to +5.25 | 4, 5, 6, 7, 8 | 0.75 to 3.50 | Gray & Brown | 15mm and 18mm | PZ ☉ ☹ 250 |
| Drivewear® | 1.5 | -6.50 to +6.00 | 4, 5, 6, 7, 8 | 0.75 to 3.50 | Drivewear® | 15mm and 18mm | DW ☉ ☹ 250 |
| DLC™ (Trivex) | 1.53 | -6.50 to +6.00 | 4, 5, 6, 7, 8, 9, 10 | 0.75 to 3.50 | Clear | 15mm and 18mm | TX ☉ ☹ 250 |
| DLC™ (Trivex) Transitions® | 1.53 | -6.50 to +6.00 | 4, 5, 6, 7, 8, 9, 10 | 0.75 to 3.50 | Gray & Brown | 15mm and 18mm | TTX ☉ ☹ 250 |
| DLC™ (Trivex) Transitions® XTRActive™ | 1.53 | -6.50 to +6.00 | 4, 5, 6, 7, 8, 9, 10 | 0.75 to 3.50 | Gray | 15mm and 18mm | TXX ☉ ☹ 250 |
| DLC™ (Trivex) Polarized | 1.53 | -6.50 to +6.00 | 4, 5, 6, 7, 8, 9, 10 | 0.75 to 3.50 | Gray & Brown | 15mm and 18mm | PTX ☉ ☹ 250 |
| 1.56 | 1.56 | -6.50 to +6.00 | 4, 5, 6, 7, 8, 9, 10 | 0.75 to 3.50 | Clear | 15mm and 18mm | 56 ☉ ☹ 250 |
| Polycarbonate | 1.586 | -6.50 to +6.00 | 4, 5, 6, 7, 8, 9, 10 | 0.75 to 3.50 | Clear | 15mm and 18mm | PC ☉ ☹ 250 |
| Polycarbonate Transitions® | 1.586 | -6.50 to +6.00 | 4, 5, 6, 7, 8, 9, 10 | 0.75 to 3.50 | Gray & Brown | 15mm and 18mm | PT ☉ ☹ 250 |

* The true base curve can be selected from a range of curves that has been defined in the Freeform Software.



ReCreating Perfect Vision®

Shamir Attitude III® Fashion

Product Availability and Markings

| Material | Index | Power Range | Nominal Base Curve* | Add | Color | Fitting Height | Invisible Marking |
|--|-------|----------------|----------------------|--------------|--------------|----------------|--|
| Polycarbonate Transitions® XTRActive™ | 1.586 | -6.50 to +6.00 | 4, 5, 6, 7, 8, 9, 10 | 0.75 to 3.50 | Gray | 15mm and 18mm | PX   250 |
| Polycarbonate Polarized | 1.586 | -6.50 to +6.00 | 4, 5, 6, 7, 8, 9, 10 | 0.75 to 3.50 | Gray & Brown | 15mm and 18mm | PPZ   250 |
| Polycarbonate Drivewear® | 1.586 | -6.50 to +6.00 | 4, 5, 6, 7, 8, 9, 10 | 0.75 to 3.50 | Drivewear® | 15mm and 18mm | PW   250 |
| Superlite™ 1.6 | 1.6 | -6.50 to +6.00 | 4, 5, 6, 7, 8, 9, 10 | 0.75 to 3.50 | Clear | 15mm and 18mm | 16   250 |
| Superlite 1.6 Tribrid™ | 1.6 | -6.50 to +6.00 | 4, 5, 6, 7, 8, 9, 10 | 0.75 to 3.50 | Clear | 15mm and 18mm | TB   250 |
| Superlite 1.6 Transitions® | 1.6 | -6.50 to +6.00 | 4, 5, 6, 7, 8, 9, 10 | 0.75 to 3.50 | Gray & Brown | 15mm and 18mm | T6   250 |
| Superlite 1.6 Transitions® XTRActive™ | 1.6 | -6.50 to +6.00 | 4, 5, 6, 7, 8, 9, 10 | 0.75 to 3.50 | Gray | 15mm and 18mm | X6   250 |
| Superlite 1.6 Polarized | 1.6 | -6.50 to +6.00 | 4, 5, 6, 7, 8, 9, 10 | 0.75 to 3.50 | Gray & Brown | 15mm and 18mm | PZ6   250 |
| Superlite 1.67 | 1.67 | -6.50 to +6.00 | 4, 5, 6, 7, 8, 9, 10 | 0.75 to 3.50 | Clear | 15mm and 18mm | 67   250 |
| Superlite 1.67 Transitions® | 1.67 | -6.50 to +6.00 | 4, 5, 6, 7, 8, 9, 10 | 0.75 to 3.50 | Gray & Brown | 15mm and 18mm | T67   250 |
| Superlite 1.67 Transitions® XTRActive™ | 1.67 | -6.50 to +6.00 | 4, 5, 6, 7, 8, 9, 10 | 0.75 to 3.50 | Gray | 15mm and 18mm | X67   250 |
| Superlite 1.67 Polarized | 1.67 | -6.50 to +6.00 | 4, 5, 6, 7, 8, 9, 10 | 0.75 to 3.50 | Gray & Brown | 15mm and 18mm | PZ67   250 |

* The true base curve can be selected from a range of curves that has been defined in the Freeform Software.



ReCreating Perfect Vision®